

Bring Sanger sequencing to the site of Ebola outbreak

Sin Hang Lee

Department of Pathology, Milford Hospital and Milford Molecular
Diagnostics Laboratory, Milford, CT, USA

Email: shlee01@snet.net

Guo Fan Hong

State Key Laboratory of Molecular Biology, Institute of Biochemistry and
Cell Biology, Shanghai Institutes for Biological Sciences, Chinese Academy
of Sciences, Shanghai China

Ebola virus disease case definition for reporting in EU

http://ecdc.europa.eu/en/healthtopics/ebola_marburg_fevers/EVDcasedefinition/Pages/default.aspx

Clinical criteria

Laboratory criteria

Any of the following:

- Detection of Ebola virus nucleic acid in a clinical specimen and confirmation by sequencing or a second assay on different genomic targets.
- Isolation of Ebola virus from a clinical specimen.

Epidemiological criteria

There are two components of an Ebola virus nucleic acid diagnostic test:

(1) Detection;


(2) Confirmation by DNA sequencing.

Is it Ebola? Malaria? Or MERS?



Perform a same-nested PCR at the site of outbreak in the shade under a tree, then the answer may become obvious.

Equipment needed to perform nested RT-PCR screen and prepare amplicons to be used as sequencing template is inexpensive



Gel electrophoresis (battery-run)



Primary PCR by pipetting



Nested PCR by microglass rod



UV-viewer (battery-run)



PCR thermal cycler




Small generator for thermal cycler

Same-nested PCR increases the detection sensitivity for RNA from Zaire Ebolavirus, Mayinga (Material generously supplied by BEI Resources www.beireseources.org)


BEI Resources Cat. No. NR-31806, Lot No. 60428456 Ebola RNA diluted in TE buffer to threshold detection level for lane 1 primary PCR, and further 10-fold serially diluted to form a decreasing ladder of concentrations for PCRs in lanes 2-7. Lane 8 was MS2 RNA positive control (band cut off)

Result: Same-nested PCR increases the NP1 and NP2 PCR detection sensitivity **by 10-fold** over a single PCR




BEI Resources Cat. No. NR-31806, Lot No. 60428456 Ebola RNA diluted in TE buffer to threshold detection level for lane 1 primary PCR, and further 10-fold serially diluted to form a decreasing ladder of concentrations for PCRs in lanes 2-7.

Result= Same-nested PCR increases the Filo A and Filo B PCR detection sensitivity **by 10-fold** over a single PCR



The non-infectious nested PCR amplicons can be safely transported to a **regional laboratory for Sanger sequencing**



Computer-generated base-calling electropherogram of a sequence of the Zai NP1/NP2-amplified 268 bp same-nested PCR product-using Zai NP2 as the sequencing primer




S/N G:1242 A:1341 T:1368 C:1576
 KB.bcp
 KB 1.2 Cap:3

DNA_2015-05-22_C01_003
 3_NP2
 KB_3130_POP7_BDTv1.mob
 Pts 1739 to 4224 Pk1 Loc:1739
 Version 5.2 Patch2 HISQV Bases: 172

Inst Model/Name 3100/LEE3130-18320-036
 May 22,2015 11:21AM, EDT
 May 22,2015 11:43AM, EDT
 Spacing:10.63 Pts/Panel1200
 Plate Name: 2015-05-22



Computer-generated base-calling electropherogram of a sequence of the FiloA/B-amplified 419 bp same-nested PCR product-using Filo A as the sequencing primer



A novel low temperature PCR makes it possible


* **Low temperature PCR** - the gateway to diagnostic Sanger sequencing for infectious diseases: **minimize *Taq* errors**

- A **HiFi** polymerase(s) for **60-120 cycles** of precision PCR amplification
- Chemical-assisted denaturing at **85°C**, instead of 94-95°C
- Chemical stabilizers for enzymes (**polymerases, ribonuclease inhibitor and reverse transcriptase**) and **dNTPs**
- Melting chemicals to **reduce mispriming**
- Wide-ranged PCR master mix **stored at 4-40°C** for clinical diagnostics
- **No pipet transferring** of post-PCR products to reduce contamination
- **No post-PCR purifications** before sequencing
- Use crude sample for primary PCR
- Adjust annealing temperature (40-50°C) for PCR stringency


* *Hong G, Lee SH, Ge S, Zhou S. A Novel Low Temperature PCR Assured High-Fidelity DNA Amplification. International Journal of Molecular Sciences. 2013; 14:12853-12862.*

sRNA from rhizobia and enterobacteria phage MS 2 RNA amplified by RT-PCR mix after 4 weeks storage at 26°C

M H2O 132bp



MS2 RT-PCR



40-cycle low temperature RT-PCR amplification of a grass carp RNA virus causing haemorrhagic disease


Target: 326 bp of viral dsRNA of grass carp reovirus.

Forward primer:


CTCAATCCTACCGGTAAGCT

Reverse primer:

ATGGTCCGCTGAAGGTTTCAT



Diagnostic base-calling electropherogram of the 326 bp of the low temperature RT-PCR amplicon of grass carp reovirus



Appenzeller T. Democratizing the DNA sequence.

Science 1990; 247:1030-2.

(PCR was invented to make templates for DNA sequencing, not for diagnosis)



Sanger sequencing of 16S rRNA gene showing mixed infections of *B. burgdorferi* and one of *B. miyamotoi* in a single patient - three characteristic double base peaks marked Z,Y,X.

Lee SH et al. Inter J Mol Sci. 2014; 15:11364-11386.





Fig. 8

```
#1 CATGCTCCACCGCTTGTGCGGGCCCCCGTCAATTCCTTTGAGTTTCACTCTTGCGAGCATACTCCCCAGGCGGCACACTTAAACACGTTAGCTTCGGTACTAACTTTT
#2 CATGCTCCACCGCTTGTGCGGGCCCCCGTCAATTCCTTTGAGTTTCACTCTTGCGAACAATACTCCCCAGGCGGCACACTTAAACACGTTAGCTTCGGTACTAACTTTT
#3 CATGCTCCACCGCTTGTGCGGGCCCCCGTCAATTCCTTTGAGTTTCACTCTTGCGAACAATACTCCCCAGGCGGCACACTTAAACACGTTAGCTTCGGTACTAACTTTT
#4 CATGCTCCACCGCTTGTGCGGGCCCCCGTCAATTCCTTTGAGTTTCACTCTTGCGAGCATACTCCCCAGGCGGCACACTTAAACACGTTAGCTTCGGTACTAACTTTT
#5 CATGCTCCACCGCTTGTGCGGGCCCCCGTCAATTCCTTTGAGTTTCACTCTTGCGAGCATACTCCCCAGGCGGCACACTTAAACACGTTAGCTTCGGTACTAACTTTT
#6 CATGCTCCACCGCTTGTGCGGGCCCCCGTCAATTCCTTTGAGTTTCACTCTTGCGAACAATACTCCCCAGGCGGCACACTTAAACACGTTAGCTTCGGTACTAACTTTT
#7 CATGCTCCACCGCTTGTGCGGGCCCCCGTCAATTCCTTTGAGTTTCACTCTTGCGAACAATACTCCCCAGGCGGCACACTTAAACACGTTAGCTTCGGTACTAACTTTT
#8 CATGCTCCACCGCTTGTGCGGGCCCCCGTCAATTCCTTTGAGTTTCACTCTTGCGAGCATACTCCCCAGGCGGCACACTTAAACACGTTAGCTTCGGTACTAACTTTT
```


Wild *Borrelia burgdorferi* isolate from a patient's blood with more than 1 copy of 16S rRNA gene

Case 2014 02 21 (template prepared by same-nested PCR after 60-cycle amplification)

To be presented at the AACC July 2015 meeting



Lower the PCR temperatures by 10°C to prepare Sanger sequencing template



You can detect infectious agents for Lyme disease, Ebola, Malaria, Cholera, Mers...at site of outbreak with

- **Extremely high specificity**
- **Extremely high sensitivity**
- **Very practical**
- **Financially sustainable**

Bring the LoTemp RT-PCR and PCR reagents in your carry-on to the field lab in West Africa (No freezing, no dry ice)

